

STATSEAL®


Patented Technology to Stop Bleeding, Absorb Exudates and Protect Wounds

What is StatSeal Powder?

StatSeal is a topical powder that forms a seal to stop the flow of blood and exudate while protecting the wound. The powder ingredients are a hydrophilic polymer and potassium ferrate.

How does StatSeal Powder work?

The mechanism of action is a simultaneous two-step action whereby the hydrophilic polymer rapidly dehydrates the blood and absorbs exudate while the potassium ferrate agglomerates blood proteins to create a seal. StatSeal does not depend on the body's natural clotting mechanism to form a seal.



Actual scanning electron microscopy image



Before StatSeal Application



With StatSeal

StatSeal stops bleeding and oozing from:

- Skin Tears
- Avulsions
- Venous Stasis Ulcers
- Pressure Ulcers
- Diabetic Ulcers
- Non-suturable Lacerations
- External Tumors
- Cancerous Lesions
- Stomas
- Other chronic or acute wounds

Benefits of StatSeal Powder:

- Creates and maintains an optimal moisture environment for wound healing
- Effective in reducing erythema and edema according to clinical studies*
- Creates a low pH environment
- Not biologically derived
- Non-systemic, non-irritant



StatSeal Powder - For External Wounds & Procedures

Manufacturer's Item

6 Sterile Bottles per Box
24 Sterile Bottles per Case

LP606
LP636



StatSeal Advanced - For Larger External Wounds

Manufacturer's Item

6 Sterile Bottles per Box
24 Sterile Bottles per Case

LP607
LP637

*Studies on file at Biolife

STATSEAL®

**Patented Technology to
Stop Bleeding,
Absorb Exudates
and Protect Wounds**


EASY TO APPLY

Clean wound. Apply StatSeal Powder when bleeding or drainage resumes.

1. Pour a thick layer of StatSeal Powder onto blood or exudate. Cover bleeding area completely. Occasionally upon application, a temporary warming sensation is felt.
2. Apply pressure for at least 2 minutes with gloved fingers or non-adherent dressing to form the protective seal. If bleeding or oozing continues, repeat application. Avoid disturbing the seal once the bleeding or oozing stops.
3. Protective seal forms. Dress wound with non-adherent dressing. The protective seal is safe to leave in place until it falls off naturally.

StatSeal Powder as applied to:


Skin Tears:



Cancerous Lesions (Squamous Cell Carcinoma)



Peristomal Ulcer:



To Remove StatSeal:
If StatSeal Powder must be removed as determined by a healthcare professional, brush away loose powder then soak wound in sterile saline until seal starts loosening from wound. Gently scrape off loosened seal with gloved finger or tongue depressor.