

STATSEAL[®]

ADVANCED

STATSEAL ADVANCED - TOPICAL POWDER

Codes: LP007 (unit)

LP607 (Box of 6 units*)

LP637 (Case of 4 boxes of 6 units*)

Country of Manufacture and Source: United States of America

Manufacturer: Biolife, L.L.C.

Description: Hemostatic topical powder for temporary external control of bleeding from sheath and catheter removals

Presentation: Application includes 1 bottle, 1 Powder Containment Device (PCD), 1 PCD foam center, and 1 product description label

Composition: Hydrophilic polymer and potassium ferrate powder.

Contraindications: None

Sterile Product, disposable and for temporary external use.**

No special storage conditions required.

Product cleared by the FDA

**The box and unit are properly marked with manufacturing lot and expiration date of the product.*

***Devices are terminally sterilized (after being placed on their respective primary packaging) with gamma radiation.*

