

STATSEAL® DISC

STATSEAL DISC – Large (~10-14 Fr)

Codes: DL016 (unit)

DL616 (Box of 6 units*)

DL626 (Case of 4 boxes of 6 units*)

Country of Manufacture and Source: United States of America

Manufacturer: Biolife, L.L.C.

Description: Hemostatic topical disc used for sealing vascular and non-vascular access sites.

Presentation: Application includes 1 StatSeal® Disc and 1 product description label.

Composition: Hydrophilic polymer and potassium ferrate powder.

Contraindications: None

Sterile Product, disposable and for temporary external use.**

No special storage conditions required.

Product cleared by the FDA

*The box and unit are properly marked with manufacturing lot and expiration date of the product.

**Devices are terminally sterilized (after being placed on their respective primary packaging) with gamma radiation.

