

STATSEAL[®]

ADVANCED

DISC

STATSEAL ADVANCED - TOPICAL POWDER

Codes: DA018 (unit)

DA618 (Box of 6 units*)

DA628 (Case of 4 boxes of 6 units*)

Country of Manufacture and Source: United States of America

Manufacturer: Biolife, L.L.C.

Description: Hemostatic topical disc for temporary external control of bleeding from sheath and catheter removal

Presentation: Application includes 1 StatSeal[®] ADVANCED Disc and 1 product description label.

Composition: Hydrophilic polymer and potassium ferrate powder.

Contraindications: None

Sterile Product, disposable and for temporary external use.**

No special storage conditions required.

Product cleared by the FDA

**The box and unit are properly marked with manufacturing lot and expiration date of the product.*

***Devices are terminally sterilized (after being placed on their respective primary packaging) with gamma radiation.*

